




Normal Pattern



Early Scleroderma Pattern



Active Scleroderma Pattern



Late Scleroderma Pattern

Images courtesy of Florian Meier and Inspectis.



Quality of Inspection

Inspectis Optical Systems
Stockholm, Sweden
www.inspect-is.com
info@inspect-is.com

Capillaroscope

Nailfold Video Capillaroscopy for Scleroderma Screening



Powered by INSPECTIS® CAP Pro Software

Delivered by world's leading expert in Capillaroscopy Microscopes and Systems!

V1, 2020. We reserve the right to modify data and specifications without prior notice.

New Generation Digital Capillaroscopy!

Cost effective

Most popular

Plug & Play

All Inclusive

Capillaroscope
Light

Capillaroscope
Self-Install

Capillaroscope
Pre-Installed


Capillaroscope
Extensive

CAP-1

CAP-2

CAP-3

CAP-4



Item No.	Item Name	CAP-1	CAP-2	CAP-3	CAP-4
IN-001	Digital Capillaroscope, USB3.0	●	●	●	●
IN-012	200x Lens	●	●	●	●
IN-024	Immersion Oil	●	●	●	●
IN-025	INSPECTIS® CAP Basics software	●			
IN-026	INSPECTIS® CAP Pro software		●	●	●
IN-033	Single-pedal Foot Switch	●			
IN-034	Triple-pedal Foot Switch		●	●	●
IN-036	Capillaroscope bench-top Holder	●	●	●	●
IN-037	Carrying Case		●	●	
IN-038	Laptop or PC			●	●
IN-039	Laptop soft case			●	
IN-040	Health care Cart				●
IN-010	15-45x Varifocal Lens with Polariser		Optional		Optional
IN-011	100x Lens		Optional		Optional
IN-035	Capillaroscope Stand with XY-stage		Optional		Optional

Non-invasive and safe method for morphologic examination and analysis of microcirculation abnormalities related to rheumatic disease, especially systemic sclerosis.

INSPECTIS® CAP Pro software offers the doctors intuitive yet powerful tools for efficient Nailfold capillaroscopy.

